

Lampade al Deuterio

Lampade da 2.000 ore

Equivalenze

Variazioni

Nuove lampade

Key: J = pre-allineate, T = con Timer, B = solo bulbo, n/a = contattateci per disponibilita'

Costruttore	Modello strumento	Codice	
		Originale	Madatec Srl
ACS	LC 750/11E		G16
ACS	LC 750/12		G16
Agilent Technologies	HP 1040, 1050VW/MW/DA		HP404
Agilent Technologies	HP 1050	79883-60002	P21
Agilent Technologies	HP 1080		DO 303
Agilent Technologies	HP 1081		DO 303
Agilent Technologies	HP 1081B		DO 303
Agilent Technologies	HP 8452A	08452-60104	DS204/05RTJ
Agilent Technologies	1100 VWD	G1314-60100	DX 224/05 J
Agilent Technologies	1100 VWD		LD-AGI-104LL
Agilent Technologies	1100 VWD		UDTA6611A
Agilent Technologies	1100 DAD		LD-AGI-105LL
Agilent Technologies	1100 DAD		UDTA3611A
Agilent Technologies	1200 VWD	G1314-60100	DX 224/05 J
Agilent Technologies	1100 DAD		LD-AGI-105LL
Alpkem	510		J53(T)
Alpkem	5010		J53(T)
Altex	100		DO 946
Altex	163		DO 903 T
Altex	164		S03
Altex	165		S03
Altex	167		S03
Aminco	DW2,		DO 902
Aminco	DW2A		DO 902
Aminco	DW2S		DO 902
Aminco	DW2C		DO 946
Aminco	DW2000		DO 946
American Hi-Tech	<i>Vedi Scientific Systems, Inc.</i>		
Amersham	<i>Vedi Pharmacia</i>		
AnalytikJena	Spekol 1200	820-60021-0	DO 249 TJ
AnalytikJena	Spekol 1300 / 1500	820-60238-0	n/a
AnalytikJena	Spekord 40/50		DO 249 TJ
AnalytikJena	Spekord S100		DO 249 TJ
AnalytikJena	Spekord 200 / 205 / 210 / 250		DO 249 TJ
AnalytikJena	S10, UV/VIS		DO 249 TJ
AnalytikJena	S5/UV		DO 249 TJ
AnalytikJena	M500		DO 249 TJ
Applied Biosystems (ABI)	120	2900-0484	J16 (T)
Applied Biosystems (ABI)	120A	2900-0484	J16 (T)
Applied Biosystems (ABI)	130	2900-0484	J16 (T)

Applied Biosystems (ABI)	130A	2900-0484	J16 (T)
Applied Biosystems (ABI)	230A	2900-0484	J16 (T)
Applied Biosystems (ABI)	270A-HAT (CZE)	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF757	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF759	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF759A	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 773	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 775	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 783	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 783 A	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 785A	2900-0484	J16 (T)
Applied Biosystems (ABI)	SF 980	2900-0484	J16 (T)
Applied Biosystems (ABI)	1000S		J16 (T)
Applied Biosystems (ABI)	all above		AB-400
Barspec	All models		DO 946
BAS	200		J53 (T)
Bausch & Lomb	<i>vedi anche Milton Roy or LDC</i>		
Bausch & Lomb	200		G05
Bausch & Lomb	700		G05
Bausch & Lomb	710		G05
Bausch & Lomb	2000		G05
Bausch & Lomb	3100		DO 903 T
Bausch & Lomb	3200		DO 903 T
Bausch & Lomb	DE 50		G05
Bausch & Lomb	Genesis 2		G31
Bausch & Lomb	Genesis 2PC		G31
Bausch & Lomb	Genesis 5		G31
Bausch & Lomb	Spectronic 21, 21D		G05
Bausch & Lomb	Spectronic 501		G05
Bausch & Lomb	Spectronic 601		G05
Bausch & Lomb	Spectronic 1001		G05
Bausch & Lomb	Spectronic 1001 plus		G05
Bausch & Lomb	Spectronic 1201		G05
Beckman	24		DO 753
Beckman	25		DO 753
Beckman	26		DO 753
Beckman	34		DO 753
Beckman	35		DO 753
Beckman	36		DO 753
Beckman	100		DO 946
Beckman	155	22 887153	R01
Beckman	163	22 947029	DO 901 T
Beckman	164		DO 753
Beckman	165	236920	DO 753
Beckman	166	239372	538706
Beckman	167	236920	DO 753
Beckman	168		538711
Beckman	5200		DO 753
Beckman	Acta		DO 753
Beckman	DB		DO 802
Beckman	DBG		DO 802
Beckman	DK		DO 802
Beckman	DU		DO 802
Beckman	DU 6	585699	DO 753
Beckman	DU 7	585699	DO 753
Beckman	DU 8	585699	DO 753

Beckman	DU 40	585699	DO 753
Beckman	DU 50	585699	DO 753
Beckman	DU 60	596791	DO 753 J6
Beckman	DU 62	596791	DO 753 J6
Beckman	DU 64	596791	DO 753 J6
Beckman	DU 65	596791	DO 753 J6
Beckman	DU 68	596791	DO 753 J6
Beckman	DU 500, 520, 530		DO 543/05 J
Beckman	DU 600	514366	514366
Beckman	DU 620	514366	514366
Beckman	DU 630	514366	514366
Beckman	DU 640	514366	514366
Beckman	DU 650	514366	514366
Beckman	DU 6000	514366	514366
Beckman	DU 7000 Series	514366	514366
Beckman	P/ACE MDQ	144-667	DS 251/05 J
Beckman	System Gold	239372	DO 753
<hr/>			
Biocad	700E		J53 (T)
Biocad	Sprint	P51085-00-5020	J53 (T)
Biocad	60		J53 (T)
<hr/>			
Biorad	206		J53 (T)
Biorad	300		J53 (T)
Biorad	1305		R19
Biorad	1306		R19
Biorad	1503		DO 946
Biorad	1790	930-6106	J53 (T)
Biorad	BioDimensions	930-6106	J53 (T)
Biorad	BioFocus 2000		J53 (T)
Biorad	BioFocus 3000		J53 (T)
Biorad	CE2000	930-6106	J53 (T)
Biorad	HPE100	930-6106	J53 (T)
Biorad	MDL 1790 / 1801		J53 (T)
Biorad	Microsampler 100		n/a
Biorad	UV-1806A	930-6132	P01
<hr/>			
Bios Analytical	Diode Array 1100	2140-0590	P17
<hr/>			
Bio-Tek Instruments	<i>vedi Kontron</i>		
<hr/>			
Biotronic	BT 3030		R19
Biotronic	BT 3035		R01
Biotronic	BT 8200		P01
Biotronic	BT 9510		J53 (T)
Biotronic	BT 9520		J53 (T)
<hr/>			
Bischoff	Lambda 1010		PR38043
	DAD3L		PR38043
<hr/>			
Buchler	all models		DO 902
<hr/>			
Buck Scientific	210 VGP	210-1052- C	n/a
<hr/>			
Camag	HPTLC Scanner	3,540,105	n/a
Camag	HPTLC Scanner II	3,540,005	n/a
<hr/>			
Camspec	302		n/a
Camspec	M302		J25
Camspec	M330		J25

Camspec	M350		J25
Camspec	M501		n/a
Carlo Erba	SP 601		DO 802
Carlo Erba	UV Vis 20		J53 (T)
Carlo Erba	UV Vis 20RS		J53 (T)
Carlo Erba	other models		DO 902
Cecil Instruments	Series 1 models		n/a
Cecil Instruments	Series 2 models		R18 (T)
Cecil Instruments	202		n/a
Cecil Instruments	272		n/a
Cecil Instruments	292		n/a
Cecil Instruments	373		n/a
Cecil Instruments	393		n/a
Cecil Instruments	505		n/a
Cecil Instruments	545		n/a
Cecil Instruments	588		n/a
Cecil Instruments	594		n/a
Cecil Instruments	595		n/a
Cecil Instruments	599		n/a
Cecil Instruments	5095		n/a
Cecil Instruments	CE 212		n/a
Cecil Instruments	CE 212A		n/a
Cecil Instruments	CE 515 Double Beam		n/a
Cecil Instruments	CE 1010		R18 (T)
Cecil Instruments	CE 1020		R18 (T)
Cecil Instruments	CE 1070		R18 (T)
Cecil Instruments	CE 1200		R18 (T)
Cecil Instruments	CE 1220		R18 (T)
Cecil Instruments	CE 2012		n/a
Cecil Instruments	CE 2021	2202-0142	R18 (T)
Cecil Instruments	CE 2112		n/a
Cecil Instruments	CE 2202		R18 (T)
Cecil Instruments	CE 2272		R18 (T)
Cecil Instruments	CE 2292		R18 (T)
Cecil Instruments	CE 4400		R18 (T)
Cecil Instruments	CE 5501		R18 (T)
Cecil Instruments	CE 5502		R18 (T)
Cecil Instruments	CE 6600		R18 (T)
Cecil Instruments	CE 6600S		R18 (T)
Cecil Instruments	CE 6602		R18 (T)
Cecil Instruments	CE 6603		R18 (T)
Cecil Instruments	CE 7200		R18 (T)
Cecil Instruments	CE 8020		R18 (T)
Cecil Instruments	CE 9020		R18 (T)
Cecil Instruments	CE 9050		R18 (T)
Cecil Instruments	CE 9080		R18 (T)
Cecil Instruments	CE 9090		R18 (T)
Ciba-Corning	UV 875		P01
Coleman	55	C055-0505	R07
Coleman	124	C055-0505	R07
Coleman	550	C055-0505	R07
Coleman	552	C055-0505	R07
Columbia Sci	SA700		Y02
Constametric	Spectromonitor III		n/a

Desaga	CD60	n/a
Dionex	<i>vedi anche Gynkotek</i>	
Dionex	AD-25	PR 39065
Dionex	CES 1	40651 J10
Dionex	CZE Capil. Elec.	J10
Dionex	DSA-1	J53 (T)
Dionex	Durrum	DO 802
Dionex	PDA-100	PR 39065
Dionex	PDA-100	LD-DIO-102LL
Dionex	VDM-1	J53 (T)
Dionex	VDM-2	40651 J10
Dionex	20001	J10
Dionex	Ultimate 3000	6,080,201 n/a
DuPont	825	DO 946
DuPont	837	DO 802
DuPont	850	DO 946
DuPont	852	R25
DuPont	860	R29
DuPont	862	R29
DuPont	870	R15
Dynamax	UV-C	J53 (T)
Dynamax	UV-D11	J53 (T)
Dynamax	UV-M	J53 (T)
EM Science	All	DO 946 J
Erma	ERC 7210	DO 946
Erma	ERC 7211	DO 946
Femto	700	DO 620
Femto	700Plus	DO 620
Femto	700S	DO 620
Femto	800XL	DO 620
GAT	LCD 500	J53 (T)
GAT	LCD 502	J53 (T)
GAT	LCD 503	J53 (T)
GAT	LCD 550	J53 (T)
GBC	All AAS Instruments	WL 23881C
GBC	LC 1205	J53 (T)
GBC	LC 1206	J53 (T)
GBC	LC 1210	J53 (T)
GBC	UV VIS 918	WL 23881C
GBC	1200	GC-400
Gilford	200	P06
Gilford	220	P06
Gilford	222	P06
Gilford	230	P06
Gilford	240	P06
Gilford	250	P06
Gilford	252	P06
Gilford	260	P06
Gilford	2000	P06
Gilford	2400	P06

Gilford	2600		P06
Gilford	Response		DO 946
Gilson	112		P06
Gilson	115	100123	DS 222 J
Gilson	116		DS 222 J
Gilson	117	100123	DS 222 J
Gilson	118	100326	DS 276 TJ
Gilson	119	100326	DS 276 TJ
Gilson	151	3/10/2026	DS 276 TJ
Gilson	152	3/10/2026	DS 276 TJ
Gilson	153	3/10/2026	DS 276 TJ
Gilson	155	100326	DS 276 TJ
Gilson	156 dual	100326	DS 276 TJ
Gilson	160 Diode Array	1G041402	DS 276 TJ
Gilson	170 Diode Array	2140-60001	DS 276 TJ
Gilson	Holochrome	30510	G09 (T)
Gilson	Spectrochrome	2900-0496	n/a
Gynkotec	UVD 160	5053.1204	DX 201 RJ
Gynkotec	UVD 320	5053.1204	DX 201 RJ
Gynkotec	UVD 170	5053.1204	DX 201 RJ
Gynkotec	UVD 340	5053.1204	DX 201 RJ
Gynkotec	UVD 170U	5053.1204	DX 201 RJ
Gynkotec	UVD 340U	5053.1204	DX 201 RJ
Gynkotec	all above		LD-GYN-100LL
Hach	DR/4000U	49513-00	DO 543/05 J
Hach	DR/5000	A23792	n/a
Hitachi	100-10		DO 946
Hitachi	100-30		DO 946
Hitachi	100-40		DO 946
Hitachi	100-50		DO 946
Hitachi	100-60		DO 946
Hitachi	101		R01
Hitachi	102		R01
Hitachi	111		R01
Hitachi	124	982-1035	DO 946
Hitachi	150-20	J 38-4017	PR38070
Hitachi	181		PR38070
Hitachi continued	200	885-1000	PR38070
Hitachi	220	885-1000	PR38070
Hitachi	220A	885-1000	PR38070
Hitachi	300	885-1000	PR38070
Hitachi	330	885-1000	PR38070
Hitachi	330LC	885-1000	PR38070
Hitachi	340	885-1000	PR38070
Hitachi	556, 557	885-1000	PR38070
Hitachi	655 A22	885-1000	PR38070
Hitachi	2000	885-1000	PR38070
Hitachi	3200	885-1000	PR38070
Hitachi	L2000 etc	885-3570	UDT0753-588
Hitachi	L-2400 UV Detector		PR38070
Hitachi	L-2420 UV-Vis		PR38070
Hitachi	L-2455 DAD		R16
Hitachi	L2500	885-3570	DO 651 MJ
Hitachi	L3000	885-3570	DO 651 MJ
Hitachi	L4000	885-3570	DO 651 MJ

Hitachi	L4200	885-3570	DO 651 MJ
Hitachi	L4225	885-3570	DO 651 MJ
Hitachi	L4250	885-3570	DO 651 MJ
Hitachi	L4500	885-3570	DO 651 MJ
Hitachi	L4720	885-3570	DO 651 MJ
Hitachi	L7400	890-2430	DO 655 J
Hitachi	L7420	890-2430	DO 655 J
Hitachi	L7450	890-2430	DO 655 J
Hitachi	L7455 DAD	890-2430	DO 655 J
Hitachi	U1000	885-1000	PR38070
Hitachi	U1100	885-1000	PR38070
Hitachi	U1500	885-1000	PR38070
Hitachi	U1800		DO 947 J
Hitachi	U2000	885-1000	PR38070
Hitachi	U2000		UDT4530932
Hitachi	U2001	885-1000	PR38070
Hitachi	U2800		DO 947 J
Hitachi	U3010	885-1000	PR38070 *
Hitachi	U3210	885-1000	PR38070 *
Hitachi	U3410	885-1000	PR38070 *
Hitachi	U3501	885-1000	PR38070 *
Hitachi	U4001	885-1000	PR38070 *
Hitachi	U4010	885-1000	PR38070 *

* alternative PR38060 - 80022397 with timer

Hi-Tech	StopFlow		n/a
IBM	LC-Detectors		DO 903 T
IBM	LC 9523 / LC 9563		DO 902
IBM	Spectrophotometer		R19
ICI	LC 1200		DO 946
ICI	LC 1205		J53 (T)
ICI	LC 1210		J53 (T)
IFCO	S500		J53 (T)
Instrument SA	JY-201		DO 802
Isco	1800		A04
Isco	V4		DO 303
Jasco	870		P01
Jasco	875		P01
Jasco	2050		R19
Jasco	Multi 320		R19
Jasco	UV 975A		P01
Jasco	Uvidec 1		R01
Jasco			
Jasco	Uvidec 2		R01
Jasco	Uvidec 4		R01
Jasco	Uvidec 4D		R07
Jasco	Uvidec 100		R01
Jasco	Uvidec 100 II		R01
Jasco	Uvidec 100 III		R19
Jasco	Uvidec 100 IV		R19
Jasco	Uvidec 100 V		R19
Jasco	Uvidec 100 VI		R19

Jasco	Uvidec 320		R01
Jasco	Uvidec 340		R01
Jasco	Uvidec 420		R01
Jasco	Uvidec 430		R19
Jasco	Uvidec 505		R01
Jasco	Uvidec 510		R01
Jasco	Uvidec 610		R19
Jasco	Uvidec 650		R19
Jasco	Uvidec 710		R01
Jasco	V530, 550		DO405J
Jobin Yvon	Monospec 103		DO 902
Joyce Loeb1	Chromoscan		DO 902
Joyce Loeb1	Scan 400		G01
Knauer	87.00	A0877	NY10256
Knauer	97.00	A0877	NY10256
Knauer	2151	A0877	NY10256
Knauer	Welchrom K-2000	A4071-1	PR18071E
Knauer	Welchrom K-2001	A4071-1	PR18071E
Knauer	Welchrom K-2500	A4071-1	PR18071E
Knauer	Welchrom K-2501	A4071-1	PR18071E
Knauer	Welchrom K-2600	A4071-1	PR18071E
Kontron	330	54-02003	DO 802
Kontron	332	91-00541A	DS 222/05 TJ
Kontron	335	91-91494	DX 222/05 TJ
Kontron	430	91-00541A	DS 222/05 TJ
Kontron	430A	91-91494	DX 222/05 TJ
Kontron	432	91-00541A	DS 222/05 TJ
Kontron	433	91-00541A	DS 222/05 TJ
Kontron	440 DAD	91-91095	DO 652/05 TJ
Kontron	535	91-91494	DX 222/05 TJ
Kontron	535 DAD	93-00636	DX 222/05 TJ
Kontron	540 DAD UV-Vis	54-02007	DX 241/05 TJ
Kontron	540+ DAD UV-Vis	54-02007	DX 241/05 TJ
Kontron	545V DAD UV-Vis	54-02007	DX 241/05 TJ
Kontron	710	54-02002	DO 904 T
Kontron	715	54-02002	DO 904 T
Kontron	720LC	54-02003	DO 905 T
Kontron	735		R19
Kontron	735LC		R19
Kontron	740LC		na
Kontron	810	54-02002	DO 904 T
Kontron	810P		G01
Kontron	820	54-02002	DO 904 T
Kontron	860	54-02002	DO 904 T
Kontron	72013C		P10
Kontron	Uvikon 430		DS 222/05 TJ
Kontron	Uvikon 432		DS 222/05 TJ
Kontron	Uvikon 600	54-02002	DO 904 T
Kontron	Uvikon 700	54-02002	DO 904 T
Kontron	Uvikon 710	54-02002	DO 904 T
Kontron	Uvikon 720	54-02002	DO 904 T
Kontron	Uvikon 720LC	54-02002	DO 904 T
Kontron	Uvikon 722		R19
Kontron	Uvikon 725		DO 946
Kontron	Uvikon 730		R19
Kontron	Uvikon 800	54-02002	DO 904 T

Kontron	Uvikon 810	54-02002	DO 904 T
Kontron	Uvikon 820	54-02002	DO 904 T
Kontron	Uvikon 860	54-02002	DO 904 T
Kontron	Uvikon 922	90-00782ST	DO 905 TJ
Kontron	Uvikon 923	90-00782ST	DO 905 TJ
Kontron	Uvikon 930	90-00782ST	DO 905 TJ
Kontron	Uvikon 931	90-00782ST	DO 905 TJ
Kontron	Uvikon 932	90-00782ST	DO 905 TJ
Kontron	Uvikon 933	90-00782ST	DO 905 TJ
Kontron	Uvikon 940	90-00782ST	DO 905 TJ
Kontron	Uvikon 941		DS 227/05
Kontron	Uvikon 942	54-02002	DO 904 T
Kontron	Uvikon 943	54-02002	DO 904 T
Kontron	Uvikon XL	96-90297	DO 907/05 TJ
Kontron	Uvikon XS, CSA	96-90297	DO 907/05 TJ
Kortek	K95	100-3609	F03
Kratos	<i>vedi Applied Biosystems</i>		
LDC-Milton Roy (TSP)	<i>vedi anche Bausch & Lomb</i>		
LDC	3000	900347	DO 903 T
LDC	3100	900347	DO 903 T
LDC	3200	900347	DO 903 T
LDC	4000	900347	DO 903 T
LDC	4100	900347	DO 903 T
LDC	4200	900347	DO 903 T
LDC	5000	900347	DO 903 T
LDC	12020		na
LDC	Spectromonitor I		na
LDC	Spectromonitor II		na
LDC	Spectromonitor III	900347	DO 903 T
LDC	Spectromonitor D	900347	DO 903 T
Lear Sigler (Teledyne Monitor Labs)			
Lear Sigler	SM 400		DO 802
Lear Sigler	SM 700		DO 802
Lear Sigler	SM 810		DO 802
Lear Sigler	SM 1000		DO 802
Lear Sigler	SM 2020		DO 802
Linear	<i>vedi anche Spectra Physics & TSP</i>		
Linear	200		J53 (T)
Linear	203		J53 (T)
Linear	204		J53 (T)
Linear	205		J53 (T)
Linear	206		J53 (T)
Linear	206PH		J53 (T)
Linear	Focus		J53 (T)
Linear	Focus 2000 High Res		J53 (T)
Linear	Focus 2000 HR		J53 (T)
Linear	UV 200		J53 (T)
Linear	All above		LI-401
LKB Biochrom	<i>vedi LKB Pharmacia</i>		
LKB Pharmacia	4050		R07
LKB Pharmacia	4054		R07
LKB Pharmacia	Ultraspec I		n/a
LKB Pharmacia	Ultraspec II		n/a

LKB Pharmacia	Ultraspec III		n/a
LKB Pharmacia	Ultraspec Plus		n/a
LKB Biochrom	<i>vedi LKB Pharmacia</i>		
LKB Pharmacia	4050		R07
LKB Pharmacia	4054		R07
LKB Pharmacia	Ultraspec I		n/a
	Ultraspec II		n/a
LKB Pharmacia	Ultraspec III		n/a
LKB Pharmacia	Ultraspec Plus		n/a
LKB Produkter	2140		n/a
LKB Produkter	2141		n/a
McPherson	632		V03
McPherson	701, 707		DO 902
Meloy Labs	SA 700	390-0039	Y02
Merck/Hitachi	<i>vedi Hitachi</i>		
Metrolab	1700		R36
Micromeritics	785		DO 902
Micromeritics	786		DO 902
Micromeritics	787		R17
Micromeritics	788		R17
Micronal	B 382		DO 506
	B582		DO 405 J
Milton Roy	<i>vedi Bausch&Lomb or LDC</i>		
Nicolet	8745		L06
Optica	10		DO 506
Oriel	26-109-002		F05
Oriel	6316		DO 946
Oriel	63162		R48
Oriel	63163		R47
Oriel	63170		J05
Perkin Elmer	AAS 1100	B0148615	DO 812 TJ
Perkin Elmer	AAS 1100(B)	B0148615	DO 812 TJ
Perkin Elmer	AAS 2100	B0148615	DO 812 TJ
Perkin Elmer	AAS 2280	00570194	R06
Perkin Elmer	AAS 2380	00570194	R06
Perkin Elmer	AAS 3030 (B)	00570194	R06
Perkin Elmer	AAS 4000	00400385	R06
Perkin Elmer	AAS 4100	B0148615	DO 812 TJ
Perkin Elmer	AAS 5000	00400385	R06
Perkin Elmer	AAS 5100	00400385	R06
Perkin Elmer	AAnalyst 600 AA	B0148615	DO 812 TJ
Perkin Elmer	AAnalyst 700 AA	B0148615	DO 812 TJ
Perkin Elmer	AAnalyst 800 AA	B0148615	DO 812 TJ
Perkin Elmer	Integral 2000		DO 946
Perkin Elmer	Integral 4000	N235 1285	DO 946
Perkin Elmer	Lambda 1	C0550505	DO 946

WWW.MADATEC.COM

Perkin Elmer	Lambda 3	C0550505	DO 946
Perkin Elmer	Lambda 3A	C0550505	DO 946
Perkin Elmer	Lambda 3B	C0550505	DO 946
Perkin Elmer	Lambda 2	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 5	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 7	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 9	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 10	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 11	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 12	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 14	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 15	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 16	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 17	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 18	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 19	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 20	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 25	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 35	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 40	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 45	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 800	B0160 917	DO 650 TJ
Perkin Elmer	Lambda 900	B0160 917	DO 650 TJ
Perkin Elmer	Lambda Array 3840	C6760011	R21
Perkin Elmer	Lambda Bio	B0160 917	DO 650 TJ
Perkin Elmer	LC-55	0271 1340	DO 946
Perkin Elmer	LC-55B	0271 1340	DO 946
Perkin Elmer	LC-65B	0271 1340	DO 946
Perkin Elmer	LC-65T	0271 1340	DO 946
Perkin Elmer	LC-75	0271 1340	DO 946
Perkin Elmer	LC-85B	0271 1340	DO 946
Perkin Elmer	LC-95	0271 1340	DO 946
Perkin Elmer	LC-135	0271 1340	DO 946
Perkin Elmer	LC-235C	0271 1340	DO 946
Perkin Elmer	LC-480 PDA	B0160917	DO 650 TJ
Perkin Elmer	M46	55-0505	R03
Perkin Elmer	M55	B0160917	DO 650 TJ
Perkin Elmer	M550	B0160917	DO 650 TJ
Perkin Elmer	M550S	B0160917	DO 650 TJ
Perkin Elmer	M551	B0160917	DO 650 TJ
Perkin Elmer	M551S	B0160917	DO 650 TJ
Perkin Elmer	M552 S/N <2030	B0160917	DO 650 TJ
Perkin Elmer	M552 S/N >2030	B0160917	DO 650 TJ
Perkin Elmer	M552A	B0160917	DO 650 TJ
Perkin Elmer	M553	B0160917	DO 650 TJ
Perkin Elmer	M559	B0160917	DO 650 TJ
Perkin Elmer	M559A	B0160917	DO 650 TJ
Perkin Elmer	Series 200 DAD, Phase I	N2922010	DS 272 TJ
Perkin Elmer	Series 200 DAD, Phase II	N2922046	na
Perkin Elmer	Series 200 UV/vis	N2900484	J16 (T)
Perkin Elmer	360	00570194	R06
Perkin Elmer	370	00570194	R06
Perkin Elmer	372	00570194	R06
Perkin Elmer	373	00570194	R06
Perkin Elmer	380	00570194	R06
Perkin Elmer	400		G20
Perkin Elmer	410		G20
Perkin Elmer	420		G20
Perkin Elmer	430		G20
Perkin Elmer	460	00570194	R06

Perkin Elmer	560	00570194	R06
Perkin Elmer	575		DO 946
Perkin Elmer	603		R06
Perkin Elmer	703		R06
Perkin Elmer	2280		R06
Perkin Elmer	2380		R06
Perkin Elmer	4000		R06
Perkin Elmer	5000		R06
PE-Hitachi	124	239-0354	DO 946 J
PE-Hitachi	139		DO 946
PE-Hitachi	200	239-0354	DO 946 J
PE-Hitachi	320	239-0354	DO 946 J
PE-Hitachi	329	239-0354	DO 946 J
PE-Hitachi	330	239-0354	DO 946 J
PE-Hitachi	340	239-0354	DO 946 J
PE-Hitachi	556	239-0354	DO 946 J
PE-Hitachi	557	239-0354	DO 946 J
Perceptive Biosystem	BIOCAD Sprint		J53 (T)
Perceptive Biosystem	BIOCAD 20		J53 (T)
Perceptive Biosystem	BIOCAD 60		J53 (T)
Pharmacia	<i>vedi LKB pharmacia</i>		
Philips	<i>vedi Unicam</i>		
Pye Unicam	<i>vedi Unicam</i>		
Saitron	Monospec B		DO 802
Scientific Systems, Inc.	500, 506 UV-Vis, 525/201		J53 (T)
Scientific Systems, Inc.	506 UV-Vis		J53 (T)
Scientific Systems, Inc.	525/201		J53 (T)
Scientific Systems, Inc.	200UV		J53 (T)
Scientific Systems, Inc.	UV 200		J53 (T)
Scinco	S-1100 Series		PR 38004A
Severn Analytical	SA 6500		J53 (T)
Severn Analytical	SA 6503		J53 (T)
Severn Analytical	SA 6504		J53 (T)
Severn Analytical	SA 6508		J53 (T)
Severn Analytical	SA 6510		J53 (T)
Shimadzu	AA-670		WL 24443A
Shimadzu	AA-670G		WL 24443A
Shimadzu	AA-6200		WL 24443A
Shimadzu	AA-6300		WL 24443A
Shimadzu	AA-6500		WL 24443A
Shimadzu	AA-6800		WL 24443A
Shimadzu	CS-920	062-65055-05	WL 24443A
Shimadzu	CS-930	062-65055-05	WL 24443A
Shimadzu	CS-9000	062-65055-05	WL 24443A
Shimadzu	D300L	062-65055-05	WL 24443A
Shimadzu	MPS-2000	062-65055-05	WL 24443A
Shimadzu	LC 3A		P06
Shimadzu	LC 4A		P06
Shimadzu	LC 6A		P06
Shimadzu	LC 2010A, 2010C	228-37401	n/a

Shimadzu	SP4		P06
Shimadzu	SPD-2A, 2AM, 2AS	670-14501-02	n/a
Shimadzu	SPD-6A	062-65056-03	DX 223
Shimadzu	SPD-6AV	062-65056-03	DX 223
Shimadzu	SPD-7A	062-65056-03	DX 223
Shimadzu	SPD-7AV	062-65056-03	DX 223
Shimadzu	SPD-10A	228-34016-02	DX 250/05 J
Shimadzu	SPD-10AV	228-34016-02	DX 250/05 J
Shimadzu	SPD-10AVP	228-34016-02	DX 250/05 J
Shimadzu	SPD-10AVVP	228-34016-02	DX 250/05 J
Shimadzu	SPD-10A, AS, AV, AVP		LD-SHI-101LL
Shimadzu	SPD-10A		UDTA6211S
Shimadzu	SPD-M10A	670-14604-00	DS 243/05 J
Shimadzu	SPD-M10AV	670-14604-00	DS 243/05 J
Shimadzu	SPD-M10AVP	228-34016-02	DX 250/05 J
Shimadzu	SPD-20A	228-34016-02	DX 250/05 J
Shimadzu	SPD-20AV	228-34016-02	DX 250/05 J
Shimadzu	SPD-M10 photo diode		n/a
Shimadzu	UV-120	200-75503-01	DO 405 J
Shimadzu	UV-120-01	200-75503-01	DO 405 J
Shimadzu	UV-120-02	200-75503-01	DO 405 J
Shimadzu	UV-140	200-75503-01	DO 405 J
Shimadzu	UV-150	200-75503-01	DO 405 J
Shimadzu	UV-160	200-75503-01	DO 405 J
Shimadzu	UV-160A	200-75503-01	DO 405 J
Shimadzu	UV-180	200-75503-01	DO 405 J
Shimadzu	UV-190	200-75503-01	DO 405 J
Shimadzu	UV-200	200-75503-01	DO 405 J
Shimadzu	UV-202	200-75503-01	DO 405 J
Shimadzu	UV-210	200-75503-01	DO 405 J
Shimadzu	UV-210A	200-75503-01	DO 405 J
Shimadzu	UV-240	200-75503-01	DO 405 J
Shimadzu	UV-250	200-75503-01	DO 405 J
Shimadzu	UV-260	200-75503-01	DO 405 J
Shimadzu	UV-265	200-75503-01	DO 405 J
Shimadzu	UV-350	200-75503-01	DO 405 J
Shimadzu	UV-365	200-75503-01	DO 405 J
Shimadzu	UV-370	200-75503-01	DO 405 J
Shimadzu	UV-625	200-75503-01	DO 405 J
Shimadzu	UV-1200 Series	200-75503-01	DO 405 J
Shimadzu	UV-min-1240	200-75503-01	DO 405 J
Shimadzu	UV-1600 Series	200-75503-01	DO 405 J
Shimadzu	UV-1700 Series	200-75503-01	DO 405 J
Shimadzu	UV-2100 Series	200-75503-01	DO 405 J
Shimadzu	UV-3000 Series	200-75503-01	DO 405 J
Shimadzu	UV-3100 Series	200-75503-01	DO 405 J
Shimadzu	All above		UDTS6601S
Shimadzu	LC2010		LD-SHI-103LL
Sigma	ZWS 11		G13
Signal	NA3 Analyzer		Y01 (T)
Speck Analytical	SA 203, 302		J53 (T)
Spectra Physics	<i>vedi anche Thermo and LDC</i>		
Spectra Physics	CE 100		J53 (T)
Spectra Physics	CE 500		J53 (T)
Spectra Physics	CE 1000		J53 (T)
Spectra Physics	CE 2000		J53 (T)

Spectra Physics	3100		DO 903 T
Spectra Physics	Focus Detector	9551-0202	J53 (T)
Spectra Physics	LC 871		n/a
Spectra Physics	Linear 200		J53 (T)
Spectra Physics	Linear 201		J53 (T)
Spectra Physics	Linear 203		J53 (T)
Spectra Physics	Linear 204		J53 (T)
Spectra Physics	Linear 205		J53 (T)
Spectra Physics	Linear 206		J53 (T)
Spectra Physics	Linear 206PH		J53 (T)
Spectra Physics	SP100		J53 (T)
Spectra Physics	SP200		J53 (T)
Spectra Physics	SP770		J16 (T)
Spectra Physics	SP970		J16 (T)
Spectra Physics	SP8200	3302-9540	DO 902 T
Spectra Physics	SP8400	3302-9540	DO 902 T
Spectra Physics	SP8430	3302-9540	DO 902 T
Spectra Physics	SP8440	3302-9540	DO 902 T
Spectra Physics	SP8450	3302-9540	DO 902 T
Spectra Physics	SP8480		DS 244 TJ
Spectra Physics	SP8480XR	2900-0487	J03
Spectra Physics	SP8490	3302-9540	DO 902 T
Spectra Physics	SP8773XR	2900-0487	J03
Spectra Physics	UV-100		J53 (T)
Spectra Physics	UV-150		J53 (T)
Spectra Physics	UV-200		J53 (T)
Spectra Physics	UV-1000		J53 (T)
Spectra Physics	UV 2000		J53 (T)
Spectronic Instruments <i>vedi Bausch&Lomb</i>			
Spectronic Unicam <i>vedi Unicam</i>			
Sun Chrom	Spectra Flow 501	200,513	DX 255/05 J
Sun Chrom	Spectra Flow 600DAD		DS 275/05 J
Suntek Science Instruments <i>vedi Scientific Systems Inc.</i>			
Thermawave	Series 3 & 5	29-013586	J90
Thermawave	Metrology System 7	29-014965	J91
Thermo Separation	Focus 2000HR	9551-0202	n/a
Thermo Separation	Focus UV3000 HR	9551-0202	n/a
Thermo Separation	Spectrachrom 100		J53 (T)
Thermo Separation	Spectrachrom 200		J53 (T)
Thermo Separation	UV 100		J53 (T)
Thermo Separation	UV 150		J53 (T)
Thermo Separation	UV 200		J53 (T)
Thermo Separation	UV 1000	9551-0023	J53 (T)
Thermo Separation	UV 2000	9551-0023	J53 (T)
Thermo Separation	UV 3000	9551-0023	J53 (T)
Thermo Separation	UV 6000LP	108052	DS 252/05 TJ
Thermo Separation	Spectro Monitor 3200		G05
Thermo Separation	Spectrofocus		J23 (T)
Thermo Spectronic <i>vedi anche Unicam</i>			
Thermo Spectronic	Aquamate		DO 544 J
Thermo Spectronic	UV 1		DO 544 J
Thermo Spectronic	UV 2		DO 544 J
Thermo Spectronic	UV 3		DO 544 J

Thermo Spectronic	UV 4		DO 544 J
Thermo Spectronic	UV1-100		DO 544 J
Thermo Spectronic	all above		J64
TJA solutions	<i>vedi Unicam</i>		
Tosoh	6040		P04
Tosoh	8010		DO 901 T
Tosoh	8011		DO 901 T
Tracor	D970A		DO 902
Transgenomic	Wave Systems 3500/4500	070024	DO651 MJ
Turner	330		G01
Turner	430		G01
Unicam	600UV		n/a
Unicam	4020	4013-164-45861	n/a
Unicam	4021	4013-164-39401	G18
Unicam	4025	4013-164-45861	n/a
Unicam	4110	4013-166-66428	A06
Unicam	4225		J53 (T)
Unicam	5600 Series		L01
Unicam	8500		L01
Unicam	8740		L01
Unicam	Crystal 200 series		R36
Unicam	Helios II	9423-UV9-002E	DO 544 J
Unicam	Helios alpha, beta	9423-UV9-002E	DO 544 J
Unicam	Helios gamma, delta	9423-UV9-002E	DO 544 J
Unicam	Helios epsilon	9423-UV9-002E	DO 544 J
Unicam	Helios (as above)		J64
Unicam	LC 3		n/a
Unicam	LC/UV		n/a
Unicam	SP 6	4013-161-05992	n/a
Unicam	SP 8	4013-163-75401	L01
Unicam	SP 30UV	4013-160-24263	B01
Unicam	SP 6350	4013-161-05992	n/a
Unicam	SP 192 (AAS)	5332-134-80131	n/a
Unicam	SP 500	4013-160-94202	A01
Unicam	SP 600	4013-161-05992	n/a
Unicam	SP 700	4013-160-24263	B01
Unicam	SP 800	4013-160-24263	B01
Unicam	SP 1700	4013-160-24263	B01
Unicam	SP 1750	4013-160-24263	B01
Unicam	SP 1800	4013-160-24263	B01
Unicam	SP 1960 (AAS)	4013-161-05992	n/a
Unicam	SP 6400	4013-161-05992	n/a
Unicam	SP 6550	4013-161-05992	n/a
Unicam	SP 7500		G05
Unicam	SP 8000	4013-160-24263	B01
Unicam	SP 8100		B01
Unicam	SP 8150	4013-163-75402	L01
Unicam	SP 8200	4013-163-75402	L01
Unicam	SP 8400	4013-163-75402	L01
Unicam	SP 8450		L06
Unicam	SP 8600		L01
Unicam	SP 8620		L07
Unicam	SP 8625		L07
Unicam	SP 8700		L06
Unicam	SP 8710		L06
Unicam	SP 8715		L06

Unicam	SP 8720		L06
Unicam	SP 8730		L06
Unicam	SP 8740		L06
Unicam	SP 8745		L06
Unicam	SP 8750		L06
Unicam	SP 8800	4013-160-24263	n/a
Unicam	SP 8820	4013-160-24263	n/a
Unicam	UV 1		DO 544 J
Unicam	UV 2		DO 544 J
Unicam	UV 3		DO 544 J
Unicam	UV 4		DO 544 J
Unicam	UV 300		DO 544 J
Unicam	UV 500		DO 544 J
Unicam	UV 550		DO 544 J
Unicam	UV (as above)		J64
Varian	AA 275	56-100218-00	R17
Varian	AA 300	56-100218-00	R17
Varian	AA 400	56-100218-00	R17
Varian	AA 475	56-100218-00	R17
Varian	AA 775	56-100218-00	R17
Varian	AA 875	56-100218-00	R17
Varian	AA 975	56-100218-00	R17
Varian	AA 1275	56-100218-00	R17
Varian	AA 1475	56-100218-00	R17
Varian	AA 3000	56-100218-00	R17
Varian	Cary 1	56-100218-00	R17
Varian	Cary 50	56-100218-00	R17
Varian	Cary 100	56-100218-00	R17
Varian	Cary 210	56-180001-00	DO 802
Varian	Cary 219 late	56-180001-00	DO 802
Varian	Cary 300		R57
Varian	Cary 400		R57
Varian	Cary 500		R57
Varian	Cary 2200		DO 802
Varian	Cary 2300		DO 802
Varian	Cary 4000		n/a
Varian	Cary 5000		n/a
Varian	Cary 6000		n/a
Varian	Prostar 310	039-161-5691	DS 228/05 TJ
Varian	Prostar 340		J53 (T)
Varian	Prostar 345		J53 (T)
Varian	Spectra AA 10		R17
Varian	Spectra AA 20		R17
Varian	Spectra AA 30		R17
Varian	Spectra AA 40		R17
Varian	SpectraAA 220 FS		R17
Varian	Spectra AA 880		R17
Varian	UV 10		DO 506
Varian	UV 50	039-161-5691	DS 228/05 TJ
Varian	UV 100	039-161-5691	DS 228/05 TJ
Varian	UV 200	039-161-5691	DS 228/05 TJ
Varian	UV 634	56-100132-00	R01
Varian	UV 635	56-100132-00	R01
Varian	UV Superscan	56-100218-00	R17
Varian	Varichrom	00-951305-00	n/a
Varian	Vista 5000	039-161-5691	DS 228/05 TJ
Varian	Vista 5500	039-161-5691	DS 228/05 TJ
Varian	75 Series (AAS)	56-100218-00	R17
Varian	2050	03-915303-99	R19

Varian	2550		P01
Varian	9050	039-161-5691	DS 228/05 TJ
Waters	480LC	WAT 099499	DO 903 T
Waters	481	WAT 099499	DO 903 T
Waters	481 LC	WAT 099499	DO 903 T
Waters	484	WAT 080357	DS 223 TJ
Waters	486	WAT 080678	PR38077
Waters	486	WAT 080678	PR39056
Waters	486		DO610/05RTJ
Waters	990	WAT 021516	DO 652/05 TJ
Waters	991	WAT 021516	DO 652/05 TJ
Waters	994 PDA	WAT 021516	DO 652/05 TJ
Waters	996 PDA	WAT 052586	DS 247/05 TJ
Waters	996, 2996 PDA		LD-WAT-104LL
Waters	2487 dual wavelength	WAT 081142	P20
Waters	2487	WAT 081142	PR39055
Waters	2487		LD-WAT-105LL
Waters	2488 dual wavelength	WAT 081142	P20
Waters	2996 PDA	WAT 052586	DS 247/05 TJ
Waters	2966	WAT 052586	PR39067
Waters	Lambda max	WAT 099499	DO 903 T
Waters	LC-Module 1		P15 (T)
Zeiss (Jena)	<i>vedi AnalytikJena</i>		
Zeiss (Oberkochen)	CLD 300		R36
Zeiss	PM 2		H 30 DS
Zeiss	PM 4		H 30 DS
Zeiss	PM 6		H 30 DS
Zeiss	PMQ II		H 30 DS
Zeiss	PMQ III		H 30 DS

Nota: Questa lista corrisponde alle nostre ultime informazioni sulla compatibilita' delle nostre lampade, siamo confidenti sulla correttezza delle informazioni.

Altresi' Madatec Srl non si assume alcuna responsabilita' su possibili eventuali errori od incongruenze presenti.