

Instrument's brand	Model	Orig. Code	Equivalent Cathodeon	Equivalent Hereaus
Agilent	See Hewlett-Packard			
Alpkem	5010		J53 (prealigned) or	
Alpkem	510		J08 bare lamp (inc T)	
Altex	155			
Altex	163		G24 (inc timer)	
Altex	164/165/167		S03 Bare Lamp no Timer	D753
Altex/Hitachi	100 (HPLC)		R07	DO946
Amersham Pharmacia	See Pharmacia			
Aminco	DW2/DW2A/DW2S		G03 F	DO902
Aminco	DW2C		R07 F	DO946
Applied Biosystems	120/120A	2900-0487	JO3 (inc timer)	DS244TJ
Applied Biosystems	130/130A/773/775	2900-0487	JO3 (inc timer)	DS244TJ
Applied Biosystems	230A			DS244TJ
Applied Biosystems	471A			DS244TJ
Applied Biosystems	757	2900-0496	J16 (inc timer)	DS244TJ
Applied Biosystems	785A/759A	2900-0484	J16 (inc timer)	DS244TJ
Applied Biosystems	783/783A	2900-0496	J16 (inc timer)	DS244TJ
Applied Biosystems	SF759	2900-0496	J16 (inc timer)	DS244TJ
Applied Biosystems	980	2900-0496	J16 (inc timer)	DS244TJ
Applied Biosystems	1000S		J16 (inc timer)	
ACS	LC750/11E		G16 A	
ACS	LC750/12		G16	
Avantes	DH2000			DS271/05J
Barspec	All Models		R07	DO946
BAS	200		J53 (prealigned) or	
BAS	200		J08 bare lamp (inc T)	
Bausch & Lomb	200		G05 A	
Bausch & Lomb	700		G05 A	
Bausch & Lomb	710		G05 A	
Bausch & Lomb	2000		G05 A	
Bausch & Lomb	DE 50		G05 A	
Bausch & Lomb	Spectronic 21		G05 A	
Bausch & Lomb	Spectronic 601		G05 A	
Bausch & Lomb	Spectronic 1001		G05 A	

Bausch & Lomb	Spectronic 1201		G05	A	
Bausch & Lomb	Genysis 5		G31	A	
Beckman	24		R03		D753
Beckman	25		R03		D753
Beckman	26		R03		D753
Beckman	34		S03		
Beckman	35		S03		
Beckman	36		S03		
Beckman	100		R07		DO946
Beckman	155		R01		DO848
Beckman	163		G24 (inc timer)		DO901T (inc timer)
Beckman	164/165		S03 bare lamp no timer		D753
Beckman	166		J37 (inc timer)	C	
Beckman	167		S03 bare lamp no timer		D753
Beckman	168		J38 (inc timer)	C	
Beckman	5200 Series		S03		
Beckman	Acta Series		R03		D753
Beckman	DB		G01		
Beckman	DBG		G01		
Beckman	DK		G01		
Beckman	DU		G01		
Beckman	DU6		S03		D753
Beckman	DU7		S03		D753
Beckman	DU8		S03		D753
Beckman	DU40		S03		D753
Beckman	DU50		S03		D753
Beckman	DU60				D753J6
Beckman	DU62				D753J6
Beckman	DU64				D753J6
Beckman	DU65				D753J6
Beckman	DU530				DO543/05J
Beckman	DU600		J31	B	DS246/05J
Beckman	DU620				DS246/05J
Beckman	DU630				DS246/05J
Beckman	DU640		J31	B	DS246/05J
Beckman	DU650		J31	B	DS246/05J
Beckman	DU800		J31	B	DS246/05J
Beckman	DU6000/7000		J31	C	DS246/05J
Beckman	System Gold		S03	F	D753
Biocad	Sprint	P5-1085-00-5020	J53 (prealigned)		DS248TJ
Biocad	700E		J53 complete assembly		DS248TJ
Biocad	See also Perceptive Biosystems				
Biorad	206		J54 inc Timer		
Biorad	300		J54 inc Timer		
Biorad	1305		R19		
Biorad	1306		R19		
Biorad	1790	930-6106	J54 inc Timer		
Biorad	BioDimensions	930-6106	J54 inc Timer		

Biorad	BioFocus 2000/3000		J53 (prealigned)	DS248TJ
Biorad	CE2000		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Biorad	HPE100		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Biorad	MDL 1790		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Biorad	Microsampler 100	930-6106	J54 inc Timer	
Biorad	UV1806A	930-6132	PO1	
Biorad	Remedi			DS246J
Bio-Tek	See Kontron			
Bios Analytical	DAD1100	2140-0590	P17	
Biotronic	BT3030		R19	
Biotronic	BT3035		R01 or R03	D753
Biotronic	BT8200		P01	
Biotronic	BT9510/9520			DS248TJ
Buchler	All models		G03	DO902
Buck Scientific	210 VGP		R44	
Camag	HPTLC Scanner		G06	
Camag	HPTLC Scanner II		G29	
Carlo Erba	SP601		G02	DO802
Carlo Erba	UV Vis 20		J53 (pre aligned)	DO902
Carlo Erba	Other Models		G03	DO902
Carlo Erba	UV2ORS		J53 (prealigned)	
Cecil Instruments	1000 Series		R18 (inc timer)	
Cecil Instruments	1070 Series			
Cecil Instruments	2000 Series		R18 (inc timer)	
Cecil Instruments	202/272/292		D02	
Cecil Instruments	3000 Series		R18 (inc timer)	
Cecil Instruments	373/393		D02	
Cecil Instruments	5000 Series		R18 (inc timer)	
Cecil Instruments	505/545/588		D02	
Cecil Instruments	594/595/599		D02	
Cecil Instruments	5095		D02	
Cecil Instruments	6000 Series		R18 (inc timer)	
Cecil Instruments	7000 Series		R18 (inc timer)	
Cecil Instruments	8000 Series		R18 (inc timer)	
Cecil Instruments	9000 Series		R18 (inc timer)	
Cecil Instruments	CE1010/CE1020/ CE1070		R18 (inc timer)	

Cecil Instruments	CE1200/CE1220		R18 (inc timer)	
Cecil Instruments	CE2012		D02	
Cecil Instruments	CE2012/CE2112		D02	
Cecil Instruments	CE2021	2202-0142	R18 (inc timer)	
Cecil Instruments	CE2112		D02	
Cecil Instruments	CE212/CE212A		D02	
Cecil Instruments	CE2202/CE2272		R18 (inc timer)	
Cecil Instruments	CE2292/CE4400		R18 (inc timer)	
Cecil Instruments	CE272		D02	
Cecil Instruments	CE292		D02	
Cecil Instruments	CE515 Double Beam		D02	
Cecil Instruments	CE5501/CE5502		R18 (inc timer)	
Cecil Instruments	CE6600/CE6600S		R18 (inc timer)	
Cecil Instruments	CE6602/CE6603		R18 (inc timer)	
Cecil Instruments	CE7200/CE8020		R18 (inc timer)	
Cecil Instruments	CE9020/CE9050		R18 (inc timer)	
Cecil Instruments	CE9080/CE9090		R18 (inc timer)	
Cecil Instruments	Series 1 Models		D02	
Cecil Instruments	Series 2 Models		R18 (inc timer)	
Ciba Corning	UV875		P01	
Coleman	55			DO946
Coleman	124		R07	DO946
Coleman	550		R07	DO946
Coleman	552		R07	DO946
Columbia Sci	SA700		Y02	
Desaga	CD60		R34	
Dionex	CES 1		J10 B	
Dionex	CZE system		J10 B	
Dionex	DSA-1		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
Dionex	Durrum		G02	DO802
Dionex	20001		J10	
Dionex	VDM-1		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
Dupont	825		R01	DO946
Dupont	837		G03	DO802
Dupont	850		R03	DO946
Dupont	852		R25	
Dupont	860		R29	

Dupont	862		R29	
Dupont	870		R15	
Dynamax	UV-C/OV-D11/UV-M		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
Erma	ERC 7210		R01 D	DO848
GBC	All AA Instruments		J07	
GBC	906		J07	
GBC	932			
GBC	1200			
GBC	1205/10		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
GBC	1206		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
GBC	1210		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
Gilford	200/220/222		G02	DO802
Gilford	230/240/250		G02	DO802
Gilford	252/260		G02	DO802
Gilford	2000/2400/2600		G02	DO802
Gilford	Response		R07	DO946
Gilson	115/116/117	10-03-26	P13 (inc timer)	DS222J
Gilson	118/119	10-03-26	P18 (inc timer)	DS276TJ
Gilson	151/152/153			DS276TJ
Gilson	160 Diode Array			DS276TJ
Gilson	170 Diode Array	2140-60001	P19	
Gilson	Holochrome	HMD 9510	G09 (inc timer)	DO902
Gilson	Spectrochrome		X02 (inc timer)	
Gynkotec	UVD320			
Gynkotec	UVD340			
Gynkotec	UVD160			
Gynkotec	UVD1705			
Gynkotec	340S DAD			
Hach	DR4000U			DO542J
Heraeus	DO222/05		P09	
Hewlett- Packard/Agilent	1040	79883- 60002		
Hewlett- Packard/Agilent	1050 VW (79853A)	79853- 60002		
Hewlett- Packard/Agilent	1050VW (79853C), 1050MW (7954A)	79883- 60002	P21	
Hewlett- Packard/Agilent	1050 DA (MWD)			
Hewlett- Packard/Agilent	1080		P05	
Hewlett- Packard/Agilent	1081		P05	
Hewlett- Packard/Agilent	1081B		P05 B	
Hewlett- Packard/Agilent	1082B		P05 B	

Hewlett-Packard/Agilent	1084		P05	
Hewlett-Packard/Agilent	1084B		P05	DO303
Hewlett-Packard/Agilent	1090 DAD Ser I & II	79883-60002	P21	
Hewlett-Packard/Agilent	1100 (G1314) VW	G1314-60100	P17	
Hewlett-Packard/Agilent	1100 (G1315A) DAD	2140-0590	P19	
Hewlett-Packard/Agilent	8450			
Hewlett-Packard/Agilent	8452A DAD			DS204/05RTJ
Hewlett-Packard/Agilent	8453		P19	
Hitachi	124		R07	DO946
Hitachi	181		R03	D753
Hitachi	100-10		R07	
Hitachi	100-40		R07	DO946
Hitachi	100-50		R07	DO946
Hitachi	100-60		R07	DO946
Hitachi	101/102/111		R01	DO848
Hitachi	150-20		R16	DO946J
Hitachi	200/220/300/330/340		R16	
Hitachi	2000/3000/4000		R16	DO946J
Hitachi	6000 / 7000		R16	
Hitachi	655A22			DO946J
Hitachi	7000 LACHROM			
Hitachi	L2500		R16	DO946J
Hitachi	L3000		R16	DO651MJ
Hitachi	L4000		R16	DO651MJ
Hitachi	L4200		R16	DO651MJ
Hitachi	L4250		R16	DO651MJ
Hitachi	L4500		R16	DO651MJ
Hitachi	L4720		R16	DO946J
Hitachi	L7400		R16	DO651MJ
Hitachi	L7450		R16	DO946J
Hitachi	U1100/U3210		R16	DO946J
Hitachi	U2000		R16	DO946J
Hitachi	U3410/U3501		R16	DO946J
Hitachi	U4001		R16	DO946J
Hitachi	UL1100		R16	
Hi-Tech	Stopflow		B03	
IBM	Detectors		G24 (inc timer)	DO902
IBM	Spectrophotometer		R19	
ICI	LC1200			DO946
IFCO	S500		J53 (prealigned) or J08 bare lamp (inc T)	DS248TJ
Instruments SA	JY-201		G02 F	DO802
Isco	1800		A01	
Jasco	870		P01	
Jasco	875		P01	
Jasco	975			
Jasco	Multi 320		R19	
Jasco	UV975A		P01	

Jasco	Uvidec 1		R01	DO848
Jasco	Uvidec 2		R01	DO848
Jasco	Uvidec 4		R01	DO848
Jasco	Uvidec 4D		R07	DO848
Jasco	Uvidec 100		R01	DO848
Jasco	Univec 100 II		R01	DO848
Jasco	Uvidec 100 III		R19	DO902
Jasco	Uvidec 100 IV		R19	DO902
Jasco	Uvidec 100 V		R19	DO902
Jasco	Uvidec 100 VI		R19	DO902
Jasco	Uvidec 320		R01	DO848
Jasco	Uvidec 340		R01	DO848
Jasco	Uvidec 420		R01	DO848
Jasco	Uvidec 430		R19	DO902
Jasco	Uvidec 505		R01	DO848
Jasco	Uvidec 510		R01	DO848
Jasco	Uvidec 610		R19	DO902
Jasco	Uvidec 650		R19	DO902
Jasco	Uvidec 710		R01	DO848
Jobin Yvon	Monospec 103		G03	DO902
Joyce LoebI	Chromoscan		G03	DO902
Joyce LoebI	Scan 400		G01	DO802
Knauer	8700/9700			
Knauer	Wellchrome			
Knauer	6100		G03	DO902
Knauer	8100		G03	DO902
Knauer	LC Photometer			
Kontron	DE430			DO902
Kontron	LC332		G34 B (inc timer)	DS222/05T
Kontron	LC335		G34 F (inc timer)	DS222/05T
Kontron	LC335	longlife		DX222/05TJ
Kontron	LC430		P09 B	DS222/05T
Kontron	LC430	longlife		DX222/05TJ
Kontron	LC432		P09 B	DS222/05T
Kontron	LC433		P09 F	DS222/05T
Kontron	440DAD			DO652/05 TJ
Kontron	535VWD	longlife		DX222/05TJ
Kontron	540DAD	longlife		DX241/05TJ
Kontron	710	D802-001678	P10	
Kontron	715	D802-001678	P10	
Kontron	LC720			DO802
Kontron	LC722			DO902
Kontron	735		R19	DO902
Kontron	735LC	9600 2746	R19	DO902
Kontron	740LC		G16 A	
Kontron	810	D802-001678	P10	DO802
Kontron	810P		G01	DO802
Kontron	820	D802-001678	P10	DO802
Kontron	860		P10	DO802
Kontron	72013C	D802-001678	P10	

Kontron	Uvikon 430		P09 B	
Kontron	Uvikon 432		P09 B	
Kontron	Uvikon 720		G01	DO905T
Kontron	Uvikon 720LC		G01	DO905T
Kontron	Uvikon 722		R19	DO902
Kontron	Uvikon 725		R07	DO946
Kontron	Uvikon 730		R19	DO902
Kontron	Uvikon 740			DO905T
Kontron	Uvikon 800		G01	DO905T
Kontron	Uvikon 810			DO905T
Kontron	Uvikon 860	D802-001678	G01	DO905T
Kontron	Uvikon 9XX			DO905TJ
Kontron	Uvikon 922			DO905TJ
Kontron	Uvikon 923			DO905TJ
Kontron	Uvikon 930		P10 (Bare Lamp)	DO905TJ
Kontron	Uvikon 931		P11	DO905TJ
Kontron	Uvikon 941			DO905TJ
Kontron	Uvikon 941 Plus			DO905TJ
Kontron	Uvikon 942			DO905TJ
Kontron	Uvikon 943			DO905TJ
Kontron	Uvikon XS			DO907/05TJ
Kontron	Uvikon XL			DO907/05TJ
Kontron	Uvikon XL Plus			DO907/05TJ
Kortek	K95		F03	
Kratos	See alao ABI			
Kratos	120/130/757/759/769	2900-0487	J03 (inc timer)	
Kratos	SF770		J03 (inc timer) / J16 (inc timer)	
Kratos	773	2900-0487	J03 (inc timer)	
Kratos	775	2900-0487	J03 (inc timer)	
Kratos	783	2900-0487	J03 (inc timer)	
Kratos	785A/980	2900-0496	J16 (inc timer)	
Kratos	873		J03 (inc timer)	
Kratos	FS970		J03 (inc timer) /	
Kratos	FS980	2900-0496	J16 (inc timer)	
LC Packings	Ultimate UV4-WL			
Lear Siegler	SM400			DO802
Lear Siegler	SM700			DO802
Lear Siegler	SM810			DO802
Lear Siegler	SM1000			DO802
Lear Siegler	SM2020			DO802
LDC	121			DO902T (inc timer)
LDC	1200			DO902T (inc timer)
LDC	1222			DO902T (inc timer)
LDC	3000		G24 (inc timer)	DO902
LDC	3100		G24 (inc timer)	DO901T (inc timer)
LDC	3200		G24 (inc timer)	DO901T (inc timer)
LDC	4000 Series		G24 (inc timer)	DO901T (inc timer)
LDC	4100		G24 (inc timer)	DO901T (inc timer)

LDC	4200		G24 (inc timer)	
LDC	12020		G06	
LDC	Spectromonitor I		G06	DO901T (inc timer)
LDC	Spectromonitor II		G06	DO901T (inc timer)
LDC	Spectromonitor III	108035	G24 (inc timer)	DO901T (inc timer)
LDC	Spectromonitor D		G24 (inc timer)	DO901T (inc timer)
Linear	See also Thermo Separation products			
Linear	UV100/200		J53 (prealigned) or J08 bare lamp (inc T)	DO248TJ
Linear	UV203		J53 (prealigned) or J08 bare lamp (inc T)	DO248TJ
Linear	UV204		J53 (prealigned) or J08 bare lamp (inc T)	DO248TJ
Linear	UV205		J53 (prealigned) or J08 bare lamp (inc T)	
Linear	UV206		J53 (prealigned) or J08 bare lamp (inc T)	DO248TJ
Linear	UV206PH		J53 (prealigned) or J08 bare lamp (inc T)	
Linear	Focus 2000 High Res		J23 (inc timer)	
Linear	Focus HR	9551-0202	J54 inc. Timer	
Linear	Focus		J53 (prealigned) or J08 bare lamp (inc T)	
LKB Biochrom	See Pharmacia			
LKB Produkter	2140/2141		FO3	
LKB Pharmacia	2151			
LKB Pharmacia	4050/4054		R07 E	DO946
McPherson	Model 701/707		G03	DO902
Meloy Labs	SA700		Y02	
Merck/Hitachi	655A/22	371140005	R16	
Merck/Hitachi	L3000	371140005	R16	
Merck/Hitachi	L4000	371140005	R16	
Merck/Hitachi	L4500	371140005	R16	
Merck/Hitachi	L4520	371140005	R16	
Merck/Hitachi	L7000/LACHROM			
Merck/Hitachi	L7400		R16	
Micromeritics	785		G17	DO902
Micromeritics	786		G17	DO902
Milton Roy	See Bausch and Lomb or LDC			
Nanometrics	216 & 9400/800		R36	

Nippon Jarrell Ash	See Thermo Jarrell Ash			
Optica	10		G03	DO902
Perkin Elmer	55E	55-0505	R03	
Perkin Elmer	341/343 Plus			St46
Perkin Elmer	141/241/241MC			St75
Perkin Elmer	272/280	0057-0194	R06	
Perkin Elmer	360/370/372/373/380	0057-0194	R06	DO848
Perkin Elmer	400		G20	
Perkin Elmer	410/420/430	0057-0194	G20	
Perkin Elmer	460	0057-0194	R06	DO848
Perkin Elmer	55/550/551/552	B016-0917	R49 (inc timer)	
Perkin Elmer	552A/553/559/559A	B016-0917	R49 (inc timer)	
Perkin Elmer	555		R30	
Perkin Elmer	575		R03	
Perkin Elmer	560	B014-8615	R06	
Perkin Elmer	603/703		R06	
Perkin Elmer	1100		P16 (inc Timer)	
Perkin Elmer	1100B/2100/4100		R06	
Perkin Elmer	2280	0047-0385	R06	DO848
Perkin Elmer	2380	0047-0385	R06	DO848
Perkin Elmer	3030	0047-0385	R06	
Perkin Elmer	4000	0047-0385	R06	
Perkin Elmer	50005100Z Leeman	0047-0385	R06	
Perkin Elmer	5000/5100/5100PC		R06	
Perkin Elmer	Aanalyst 600 & 700		R06	
Perkin Elmer	Aanalyst 800		R06	
Perkin Elmer	Integral 2000		R07	DO946
Perkin Elmer	Integral 4000		R07	DO946
Perkin Elmer	Lambda 1/Lambda 3	C055-0505	R03	
Perkin Elmer	Lambda 2/2S/10/11	B016-0917	R49 (inc timer)	DO650TJ
Perkin Elmer	Lambda 3A/3B/6	C055-0505	R03	
Perkin Elmer	Lambda 3B		R12	
Perkin Elmer	Lambda 4/14/18		R49 (inc timer)	
Perkin Elmer	Lambda 5/7/9	B016-0917	R49 (inc timer)	DO650TJ
Perkin Elmer	Lambda 12/16/20/40	B016-0917	R49 (inc timer)	DO650TJ
Perkin Elmer	Lambda 15/17/19	B016-0917	R49 (inc timer)	DO650TJ
Perkin Elmer	Lambda 800/900	B016-0917	R49 (inc timer)	DO650TJ
Perkin Elmer	Lambda Bio	B016-0917	R49 (inc timer)	
Perkin Elmer	Lambda Array 3840		R21	
Perkin Elmer	LC 55	087-1171	R03	DO946
Perkin Elmer	LC 55B	087-1171	R03	DO946
Perkin Elmer	LC 65T	087-1171	R07	DO946
Perkin Elmer	LC 75	0271-1340/0271-1706	R07	DO946
Perkin Elmer	LC 85/85B	0271-1340/0271-1706	R07	DO946
Perkin Elmer	LC 95	0271-1709	R07	DO946
Perkin Elmer	LC 135		R07	DO946
Perkin Elmer	LC 235		R42	DO946
Perkin Elmer	LC 290			
Perkin Elmer	LC480			DO650TJ
Perkin Elmer	LC1355	N235-1125	R42	
Perkin Elmer	M46	55-0505	R03	
Perkin Elmer	M55	55-0505	R03	DO650TJ

Perkin Elmer	M57	55-0505	R03	
Perkin Elmer	M550	102-014	R13	DO650TJ
Perkin Elmer	M550S	092-941	R30	
Perkin Elmer	M551	102-014	R13	DO650TJ
Perkin Elmer	M551S	092-941	R30	
Perkin Elmer	M552 from S/N 2030	102-014	R13	DO650TJ
Perkin Elmer	M552 to S/N 2030	092-941	R30	
Perkin Elmer	M554	092-941	R30	
Perkin Elmer	M555	102-014	R13	DO650TJ
Perkin Elmer	M559			DO650TJ
Perkin Elmer	M559A			DO650TJ
Perkin Elmer	Series 200 DAD			DS272TJ
Perkin Elmer	Zeeman 5000		R06	
Perkin Elmer	Zeeman 5100PC		R06	
Perkin Elmer/ABI	See Applied Biosystems			
PE/Hitachi	LC90			
PE/Hitachi	124	239-0354	R16	
PE/Hitachi	139	98201035	R07	
PE/Hitachi	200	239-0354	R16	DO946J
PE/Hitachi	320	239-0354	R16	DO946J
PE/Hitachi	330	239-0354	R16	DO946J
PE/Hitachi	340	239-0354	R16	DO946J
PE/Hitachi	556	239-0354	R16	DO946J
PE/Hitachi	557	239-0354	R16	DO946J
Perseptive Biosystems	BioCAD/SPRINT system for Perfusion Chromatography		J53 (prealigned) or	
			J08 bare lamp (inc T)	
Perseptive Biosystems	Biocad 20		J53 (prealigned) or	
			J08 bare lamp (inc T)	
Pharmacia	Ultrospec II from s/n6598		R07 E	
Pharmacia	Ultrospec III		R07 E	
Pharmacia	TDS from s/n 206		R07 E	
Pharmacia	Ultrospec Plus from s/n176		R07 E	
Philips	See Unicam			
Pye	See Unicam			
Saitron	Monospec B		G01	DO802
Scientific Systems Inc (SSI)	200UV		J53 (prealigned)	DS248TJ
Scientific Systems Inc (SSI)	500		J53 (prealigned)	DS248TJ
Scientific Systems Inc (SSI)	506 UV/Vis		J53 (prealigned)	DS248TJ
Scientific Systems Inc (SSI)	525.201		J53 (prealigned)	DS248TJ
Scientific Systems Inc (SSI)	UV200		J53 (prealigned)	DS248TJ
Secomam	ASD5000		J81 (prealigned)	
Secomam	Dathelie		J81 (prealigned)	

Secomam	Pastel-UV		J81 (prealigned)	
Severn Analytical	6510		J53 (prealigned)	DS248TJ
Severn Analytical	6500		J53 (prealigned)	DS248TJ
Severn Analytical	6503		J53 (prealigned)	DS248TJ
Severn Analytical	6504		J53 (prealigned)	DS248TJ
Severn Analytical	6508 (Rapidscan)		J53 (prealigned)	DS248TJ
Shimadzu	625/630/640/670/680			DO405J
Shimadzu	6500AA		P06	
Shimadzu	CS930			DO405J
Shimadzu	CS9000			DO405J
Shimadzu	D300L		G03	DO902
Shimadzu	LC3A		P06	
Shimadzu	LC4A		P06	
Shimadzu	LC 6A		P06	
Shimadzu	SPD 2A/AV		P06	
Shimadzu	SPD 6A		P06	
Shimadzu	SPD 6AV		P06	
Shimadzu	SPD 7AV			
Shimadzu	SP 4		P06	
Shimadzu	SPD10A/SPD10AS	C228-34016-02	J74	DS250/05J
Shimadzu	SPD10AV/SPD10AVP	C228-34016-02	J74	DS250/05J
Shimadzu	SPDM10AVP		J74	DS250/05J
Shimadzu	SPD10A	longlife		DX250/05T
Shimadzu	UV120/140/150/160/180/190/200/202/210A			DO405J
Shimadzu	UV160A			DO405J
Shimadzu	UV240	200-75503-01		DO405J
Shimadzu	UV250			DO405J
Shimadzu	UV260			DO405J
Shimadzu	UV365			DO405J
Shimadzu	UV1200			DO405J
Shimadzu	UV1600			DO405J
Shimadzu	UV2100			DO405J
Shimadzu	UV3000/3100			DO405J
Sigma	ZWS 11		G13 (inc timer)	
Signal	NA3 Analyzer		Y01 (inc timer)	
Speck Analytical	SA203		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Speck Analytical	SA 302		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Spectra Physics	100/200		J53 (prealigned) or	DS248TJ
			J08 bare lamp (inc T)	
Spectra Physics	3100		G24 (inc timer)	
Spectra Physics	Focus Detector		J53 (prealigned) or	DS248TJ

		J08 bare lamp (inc T)	
Spectra Physics	LC871	A04	
Spectra Physics	SP8200	G03 F	
Spectra Physics	SP8400	G13 (inc timer) C	DO902T (inc timer)
Spectra Physics	SP8480	G13 (inc timer) C	DO902T (inc timer)
Spectra Physics	SP8480XR	J03 (inc timer)	
Spectra Physics	SP8430	G13 (inc timer) C	DO902T (inc timer)
Spectra Physics	SP8440	G13 (inc timer) C	DO902T (inc timer)
Spectra Physics	SP8450	G13 (inc timer) C	DO902T (inc timer)
Spectra Physics	SP8490	G13 (inc timer) C	
Spectra Physics	SP8773XR	J03 (inc timer)	
Spectra Physics	Spectrachrom 100	J53 (prealigned) or	DO248TJ
		J08 bare lamp (inc T)	
Spectra Physics	Spectrachrom 200	J53 (prealigned) or	DO248TJ
		J08 bare lamp (inc T)	
Spectra Physics	UV 100, 150, 200	J53 (prealigned) or	DO248TJ
		J08 bare lamp (inc T)	
Spectra Physics	UV1000/2000/3000	J53 (prealigned) or	DO248TJ
		J08 bare lamp (inc T)	
Spectra Physics	UV6000		DS252/05TJ
Spectronic	See Bausch and Lomb		
Spectronic	Genesys 2	G31 A	
Spectronic	Genesys 2PC	G31 A	
Spectronic	Genesys 5	G31 A	
Spectronic	21D	G05 A	
Spectronic Unicam	See Unicam		
Sunchrom	SF600 DAD		DS275/05J
Suntek Science Instruments	See Scientific Systems (SSI)		
Thermo Quest	See Thermo Separation and Spectra Physics/Linear		
Thermo Finnigan	Surveyor	J53 (prealigned)	
Thermo Separation	3200	J53 (prealigned)	
Thermo Separation	SC100/200	J53 (prealigned) or	DO248TJ
Thermo Separation		J08 bare lamp (inc T)	

Thermo Separation	Spectra 100		J53 (prealigned) or	DO248TJ
Thermo Separation			J08 bare lamp (inc T)	
Thermo Separation	UV100/150/200		J53 (prealigned) or	DO248TJ
Thermo Separation			J08 bare lamp (inc T)	
Thermo Separation	UV1000/2000/3000		J53 (prealigned) or	DO248TJ
Thermo Separation			J08 bare lamp (inc T)	
Thermo Separation	Focus 2000 HR	A4903-009	J54 inc. Timer	
Thermo Separation	UV3000 HR	A4903-010	J54 inc. Timer	
Thermo Separation	UV6000	24198D	J35 inc Timer	DS252/05TJ
Thermo Separation	Spectrofocus		J23 bare lamp inc T	
Thermo Spectronic	See Unicam			
Thermo Spectronic	Aquamate		J64	
Thermo Spectronic	UV1/UV2/UV3/UV4		J64	
Thermo Spectronic	UV1-100		J64	
TJA Solutions	See Unicam			
Tosoh	6040		P04	
Tosoh	8010		P02 inc. Timer	
Tosoh	8011		P02 inc. Timer	
Tracor	D970A			DO902
Turner	330		G01	DO902
Turner	430		G01	DO902
Unicam	200			
Unicam	600UV	4013-161-05993	D01	
Unicam	4020	4013-164-45861	A04	
Unicam	4021	4013-164-39401	G18	
Unicam	4025	4013-164-45861	A04	
Unicam	4110	4013-166-66428	A06	
Unicam	4225		J08 (inc timer)	DS248TJ
Unicam	8500		L01	
Unicam	8740		L01	
Unicam	5600 Series		L01 B	
Unicam	Helios II		J64	DO544J
Unicam	Helios Alpha/Beta		J64	DO544J
Unicam	Helios Delta/Gamma		J64	DO544J
Unicam	Helios Epsilon		J64	DO544J
Unicam	LC/UV	4013-164-45861	A04	
Unicam	LC3	4013-164-45861	A04	
Unicam	SP6	4013-161-05992	D01	

Unicam	SP6-350	4013-161-05992	D01	
Unicam	SP8		L01	
Unicam	SP30 UV	4013-160-24263	B01	
Unicam	SP192	5332-134-80131	D01	
Unicam	SP500	4013-160-94202	A01	
Unicam	SP600	4013-161-05992	D01	
Unicam	SP700	4013-160-24263	B01	
Unicam	SP800	4013-160-24263	B01	
Unicam	SP1700	4013-160-24263	B01	
Unicam	SP1750	4013-160-24263	B01	
Unicam	SP1800	4013-160-24263	B01	
Unicam	SP1960	4013-161-05992	D01	
Unicam	SP6400	4013-161-05992	D01	
Unicam	SP6550		D01	
Unicam	SP7500	4013-160-24263	G01	
Unicam	SP8000		B01	
Unicam	SP8100		B01 or	
Unicam	SP8100	4013-163-75402	L01	F
Unicam	SP8150	4013-163-75402	L01	
Unicam	SP8200	4013-163-75402	L01	
Unicam	SP8400		L01	
Unicam	SP8450		L06	B
Unicam	SP8600		L01	B
Unicam	SP8620		L07	B
Unicam	SP8625		L07	B
Unicam	SP8700		L06	B
Unicam	SP8710/8715		L06	B
Unicam	SP8720/8730/8740		L06	B
Unicam	SP8745		L06	B
Unicam	SP8750	4013-160-24263	L06	B
Unicam	SP8800	4013-160-24263	B01	B
Unicam	SP8820		B01	B
Unicam	UV2		J64	DO544J
Unicam	UV300/500/550		J64	DO544J
Uvikon	See Kontron			
Varian	UV 100/200	56-100218-00		
Varian	75 Series	03-915303-99	R17	
Varian	2050		R19	

Varian	2550		P01	
Varian	5500	56-100218-00		
Varian	AA 275	56-100218-00	R17	
Varian	AA 400	56-100218-00	R17	
Varian	AA 475	56-100218-00	R17	
Varian	AA 775	56-100218-00	R17	
Varian	AA 875	56-100218-00	R17	
Varian	AA 975	56-100218-00	R17	
Varian	AA1275	56-100218-00	R17	
Varian	AA1475	56-100218-00	R17	
Varian	AA 3000	56-180001-00	R17	
Varian	Cary 100		R17	
Varian	Cary 210	56-180001-00	G01	
Varian	Cary 219 (later)	56-180001-00	G01	
Varian	Cary 500		R17	
Varian	Cary 2200	56-180001-00	G01	
Varian	Cary 2300	56-180001-00	G01	
Varian	ProStar 340/345		J53 (prealigned)	DS248TJ
Varian	Spectra A10	56-180001-00	R17	
Varian	Spectra A20	56-180001-00	R17	
Varian	Spectra A30	56-180001-00	R17	
Varian	Spectra A40		R17	
Varian	UV10	00-951305-00	G03	
Varian	UV100			
Varian	UV50	56-100132-00	G14	
Varian	UV634	56-100132-00	R01	
Varian	UV635	56-100218-00	R01	
Varian	UV Superscan	00-951305-00	R17	
Varian	Varichrom		G14	
Varian	9050			
Waters	450			
Waters	480/480LC		G24 (inc timer)	DO903T (inc timer)
Waters	481		G24 (inc timer)	DO903T (inc timer)
Waters	481LC		G24 (inc timer)	DO903T (inc timer)

Waters	Lambda Max	80357	G24 (inc timer)	DO903T (inc timer)
Waters	484	80678 Except Vr 4.0 or above software)	P15 (inc timer)	DS223TJ
Waters	486		P15 (inc timer)	DO610/O5RTJ
Waters	455			
Waters	487 (aka 2487)		P20	
Waters	990			DO652/05TJ
Waters	991			DO652/05TJ
Waters	994 DAD			DO652/05TJ
Waters	996		J66	DS247/05TJ
Waters	UV2487/Alliance	WAT 081142	P20	DS225/05J
Waters	Alliance	Check model number of detector in use		
Waters	LCI		P15 (inc timer)	
Zeiss	CID 300		R36	
Zeiss	PM6			H30DS
Zeiss	PMQ11/PMQ111			H30DS
Zeiss	M4Q3			H30DS
Zeiss	Spekol			DO249TJ
Zeiss	Specord			DO249TJ
Zeiss	S10			DO249TJ
Zeiss	UV/Vis			DO249TJ
Zeiss	S5UV			DO249TJ
Zeiss	M500			DO249TJ