

Xenon and Mercury

Costruttore	Modello dello strumento	Codice
Agilent	1046A	2140-0549
Aminco-Bowman	Fluorescence (most types)	XBO 150W/S
Aminco-Bowman	SLM4800 Ushio Lamp	UXL 450S-O
Baird Atomic	Fluorimipoint/FP100, RC200 FP; Fluoricord/FC100, RC200 FP; Fluoricord/FC100, RC200 FC; Fluorispec; SF1: SF100E; SFR100; Nova	XBO 150W/S
Cary	60 CD ORD	XBO 450W/4
Ciba Corning (see Jasco)		
Farrand	FOCI Mk 1, Mk 2, System 3	XBO 150W/S
Hitachi	F2000; F4010	XBO 150W/1
Hitachi	F1000; F1050; F1080 (Standard)	XBO 150W/1
Hitachi	(For enhanced UV)	XBO 150W/4
Jasco	SP820; SP821; SP210; 820FP; 821FP, 821, 920	XBO 150W/1
Jasco	920 Ushio Lamp	UXL-150S
Kontron	SFM 25	XBO 150W/1
LDC (TSP)	FM4100	XBO 150W/1
Leitz Microscopes	50 Watt Versions	HBO 50W AC/L1
Leitz Microscopes	100 Watt Versions	HBO 100W/2
Merck Hitachi	Various. Ushio Lamp	UXL-152H
Merck Hitachi	F1000; F1050; F1080 Std	XBO 150W/1
Merck Hitachi	(Various Lamp Options) Enhanced UV	XBO 150W/4
Nikon	Diaphot Fluorescence Microscope	HBO 100W/2
Nikon	Optiphot 2; uses various lamps	HBO 100W/2
		XBO 100W/Ofr
Olympus	IMT-2 Microscope	HBO 100W/2
Linear	LC304 (Lamp on Bracket)	9551-0144
Perkin-Elmer	203; 204; MPF 2A; MPF3; MPF4; 650;hplc 200/200A	XBO 150W/1
Philips	PU4027	XBO 150W/1
Photon Technology	DeltaScan	XBO 75W/2
Shimadzu	RF 540; RF 535: RF 551	XBO 150W/1-SM
Shimadzu	10 AXL. Ushio Lamp	UXL-150S
Spectra-Physics (TSP)	FL2000; LC304, Ultrafluor, FL3000	9551-0144
Spex	HPI HPLC Monitor	XBO 150W/S
Varian	9070,00	9551-0144
Vickers	55 Projection Microscope: M55	XBO 150W/1
Waters	470, 474. Ushio Lamp	UXL-150MO