

Costruttore e modello strumento	Codice originale	Madatec Part Number	Note
ABI	2450-0213		
AGFA	F7.9924.4236	CT100T03	
Agilent			see Hewlett Packard
Bayer	071-0392-01	CT100T03	
Bayer Immuno	156-2602-01	CT065T01	
Beckman	239370		-
Beckman 166	M47	CT020T01	12v/20W bare lamp only
Beckman DU7	873947		12v/50W
Beckman DU70		CT020T01	12v/20W bare lamp only
Beckman DU600	945672	CT020T01	
Beckman Coulter		See Beckman	
Biochrom		See Pharmacia	
BLV	115320	CT020T01	12v/20W
BLV	115520	CT050T01	12v/50W
BLV	115530	CT050A01	12v/50W
BLV	116910	CT010T03	6v/10W
Camag	TLC Scanner 3		
Cathodeon	N/A	CT035A01	12v/35W
Cathodeon	A1/223 or EHJ	CT250T01	24v/250W
Cathodeon	A1/239/EVD	CT400T01	36v/400W
Cathodeon	AF5/4	TBA	TBA
Cathodeon	EXL	CT030A02	6.6A/30W
Cathodeon	J1/39	CT200T01	6.6A/200W
Cathodeon	J1/57	CT045A01	6.6A/45W
Cathodeon	J1/58	CT100T04	6.6A/100W
Cathodeon	J1/59	CT036A01	6A/36W
Cathodeon	J1/65	CT200T02	8.33A/200W
Cathodeon	J1/66	CT100T01	6.6A/100W
Cathodeon	J1/74	CT030A01	6.6A/30W
Cathodeon	M28 or EVA	CT100T02	12v/100W
Cathodeon	M29	CT010T02	6v/10W
Cathodeon	M32	CT050T01	12v/50W
Cathodeon	M42	CT010T03	6v/10W
Cathodeon	M47	CT020T01	12v/20W
Cathodeon	M74	CT050A01	12v/50W
Cathodeon	M75	CT035A01	12v/35W
Cathodeon	M86	CT020A01	12v/20W
Cecil 200 Series		CT035A03	12v/35W
Cecil 272		CT035A03	12v/35W
Cecil 393		CT035A03	12v/35W
Cecil CE 1010/CE 1020			Tungsten Halide Lamp
Cecil 2000 Series		CT035A03	12v/35W
Cecil CE 2041			
Cecil 2292		CT035A03	12v/35W
Cecil 3000 Series		CT035A03	12v/35W
Cecil TH Lamp	2202-0140	CT035A03	12v/35W
Cecil TH Lamp	2303-0140	CT035A03	12v/35W

Costruttore e modello strumento	Codice originale	Madatec Part Number	Note
GE Q10T2.5	M29	CT010T02	6v/10W
GE Q50T3	M32	CT050T01	12v/50W
GE Q20T2.5	M47	CT020T01	12v/20W
GE Q50T03	M74	CT050A01	12v/50W
GE Q35T3	M75	CT035A01	12v/35W
GE Lighting	M29	CT010T02	6v/10W
GE Lighting	M32	CT050T01	12v/50W
GE Lighting	M42	CT010T03	6v/10W
GE Lighting	M47	CT020T01	12v/20W
GE Lighting	M74	CT050A01	12v/50W
GE Lighting	M75	CT035A03	12v/35W
Gilson 118	100325		
Hewlett Packard 8450A	08450-60105		
Hewlett Packard 8453	G1103-60001		
Hewlett Packard 8452A NIR	08452-60121		
Hewlett Packard	G1103-600001		
Hewlett Packard	HP1100 G1315A DAD		
Hitachi	J851170		10V/4A Tungsten Lamp
Hitachi	84-0300-69		10V/4A Tungsten Lamp
Hitachi	239-3018		10V/4A Tungsten Lamp
Hitachi L4200	885-1200		Tungsten lamp Lamp
Jasco	V530/550/560/570		
Jasco	5330-0087		
Jasco	7800		
Jenway	605-064		
Jenway 6501 Colorimeter			
Jenway 6500			
Kernchen 7387/7388	M30/ESB		6v/20W
Knauer	A4072		
Milton Roy 2000	333323		24v/24W pre-aligned
Narva (aka GLE)	YT55134B		12v/20W pre-aligned
Orbitec	H64425/M47	CT020T01	12v/20W
Orbitec	130108/M47	CT020T01	12v/20W
Osram 64663 HLX	A1/239/EVD	CT400T01	36v/400W
Osram 64655 HLX	A1/223/EHJ	CT250T01	24v/250W
Osram 64386	J1/39	CT200T01	6.6A/200W
Osram 64311	J1/59	CT036A01	6A/36W
Osram 64321	J1/57	CT045A01	6.6A/45W
Osram 64346	J1/58	CT100T04	6.6A/100W
Osram 64347	J1/66	CT100T01	6.6A/100W
Osram 64387	J1/65	CT200T02	8.33A/200W
Osram 64225	M29	CT010T02	6v/10W

Costruttore e modello di strumento	Codice originale	Madatec Part Number	Note
Osram 64425 (aka Halostar) Osram	M47 64258	CT020T01	12v/20W
PE/Hitachi 557	557-0240		Pre-aligned
Perkin Elmer Perkin Elmer Perkin Elmer 55, 55SA, 550, 551, 552, 553, 559 & 559A	3009-1906 B007-2904 B009-1906	CT030T01	
Perkin Elmer Lambda 2, 5, 7, 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 40, 900 & Bio	B011-4620	CT030T02	Pre-aligned
Perkin Elmer Lambda 1 & 3 Perkin Elmer 550A	C055-0500 C055-0500		Tungsten Halide lamp
Perkin Elmer Lambda 4 & 6 Perkin Elmer Perkin Elmer DAD 200	C688-0054 M551S N2922011	PED200	
Pharmacia Pharmacia Pharmacia 2000 series Pharmacia 3000 series Pharmacia Ultrospec 4050	41-8200-30 80-2106-16 M86 M86	CT020A01 CT020A01 CT020A01 CT020A01	Bare lamp Bare lamp Bare lamp only Bare lamp only
Philips Philips Philips Philips Philips Philips 7387/7388	5972/M43 13078 13102 13103 LN 29 840 M30/ESB	CT020A01 CT050A01 CT035A03	6v/10W 12v/20W 12v/50W 12v/35W 28v/85W Linear 6v/20W
Radium Radium	RJ1006/M42 RJ2012/M86	CT020A1	6v/10W 12v/20W
Roper Scientific	8401-107-17	CT250T01	
Shimadzu Shimadzu SPD10AV	062-65004-06 670-14602 or 6598379512		12v/20W Pre-aligned
Spectra Physics		See TSP	
Spectronic Instruments 20D+			
Thermo Separation Products		See TSP	
Thermo Spectronic Thermo Spectronic Thermo Spectronic	9423 UV9 0004E Aquamate UV1-100	CT020A01 CT020A01 CT020A01	Bare lamp only Bare lamp only Bare lamp only
Thermo Spectronic		See also products listed under Unicam	

Costruttore e Modello di strumento	Codice originale	Madatec Srl Part Number	Note
TJA Solutions		See Unicam	
TSP	108053	CT020T01	12v/20W
TSP UV3000	9551-0111		
Unicam	4013-232-02001	CT035A01	Bare lamp
Unicam	4013-232-02021	CT020A01	Bare lamp
Unicam UV2	9423-UV9-0001E	CT035A01	Bare lamp
Unicam 5625	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 5675 (Vis)	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 8600 Series	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 8625	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 8675 (Vis)	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 8700 Series	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam 8800 Series	9423-185-03021	CT030T03	Pre-aligned lamp
Unicam SP6	9435-240-49431		Prealigned
Unicam SP8	2413-533-90005		12v/50W
Unicam SP500/800	4013-170-00011		
Unicam SP1700/1800	4013-170-00011		
Unicam SP7500	4013-229-57161		24v/24W Pre-aligned
Unicam SP8000	4013-170-00011		
Unicam Helios Epsilon	9423-UVE-3353E	CT010T03	
Unicam Helios Gamma	9423-UV9-0004E	CT020A01	Bare lamp only
Unicam Helios Delta	9423-UV9-0004E	CT020A01	Bare lamp only
Varian	56-101332-00	CT055A02	12v/55W OEM lamp