

<p>Tube, pyrocoated, with rings</p> <p>recommended standard tube</p> <p>SYS 2000; SYS 3000</p>	99-0059-00	<p><i>GBC-Pyro-Rohr</i> Einspritzloch Injection hole</p> <p>Pyro-Schicht pyro coating</p> <p>Ringe rings</p>	56GB700
<p>Tube, coated, halfrings</p> <p>low cost, temperature program may be altered</p> <p>SYS 2000; SYS 3000</p>	Not available		56GB701
<p>Tube, pyrocoated, Omega platform</p> <p>replaces mono-snap platform</p>	Not available		56GB710
<p>Tube, uncoated, rings</p> <p>standard tube, for Pb, Cd,... NOT RECOMMENDED!</p> <p>SYS 2000; SYS 3000</p>	45-0005-00		56GB720
<p>Tube, pyrocoated, rings</p> <p>LONGLIFE, thicker coating, suitable for over night runs</p>	Not available		56GB725
<p>Platform, L'vov type</p> <p>Pyro graphite, insert in standard tube</p>	99-0060-00		56GB750
<p>Tube, for Ultra-Z</p> <p>pyro coated</p>	45-0012-00		56GB1UZ
<p>Tube, for Ultra-Z, holes and rings</p> <p>pyro coated</p>	Not available		56GB2UZ
<p>Electrode / Contact</p> <p>set of 2 pieces</p>	99-0061-00		56GB811
<p>Shroud</p> <p>price per unit</p>	45-0004-00		56GB902