









Tube, pyrocoated, with segment	180-7444-ELC	 D7,8/4,7x30mm	56Hi001
Tube, uncoated, with segment	180-7400	 D7,8/4,7x30mm	56Hi002
Tube, pyrocoated, with segment for extended injection volume	190-6003		56Hi003
Platform, L'vov type NOT RECOMMENDED!	180-7404	 D18x3,8x1,2 mm	56Hi004
Tube, pyrocoated, for forked platform	190-6007 (ANO 6007)	 D8/5,0x28mm	56Hi005
Forked Platform Pyro graphite, insert in 56Hi005, NOT RECOMMENDED	190-6008 (ANO 6008)		56Hi006
Tube, coated, with forked platform	190-0028 (ANO 0028)		56Hi007
Electrode / Contact set of 2 pieces	180-7401	 D20/4x21,5 mm	56Hi008