

The Minidiff Plus Diffuse Reflectance Accessory For FT-IR Spectroscopy



The Minidiff Plus diffuse reflectance accessory is designed for routine sampling. High performance solids analysis is made simple through the optimised zero alignment optics and Benchmark™ baseplate mounting system. Sampling versatility is maintained through the choice of either standard solids cup or Diabrase™ abrasive sampling pads. Up to 3 samples can be loaded at the same time, while the sample introduction system reduces the risk of sample spillage.

Features

Zero Alignment

- Completely pre-aligned with simple sample introduction/removal.

Optimised Throughout

- Designed for a high throughput on all FT-IR spectrometers.

Maximum Versatility

- Sample cups for powders and 'Diabrase' sampler for non particulate solid samples.

Budget Cost

- Cost engineered to meet the budget of all laboratories.

Ordering Information

4510 Minidiff Plus diffuse reflectance accessory.

Includes: Optics system with Benchmark™ instrument baseplate, abrasive sample holders (2), sample cup holders (2), Diabrase™ abrasive pads (20), KBr powder (50g) and pestle and mortar.

Please specify instrument make and model.

4505 Spare sample cup holders with 3 cups (2)

4506 Spare packet of Diabrase™ sampling pads (100)

4508 Spare abrasive sample holders with 3 pads (2)

3610 KBr powder (50g)

3600 Pestle and mortar

THE SOLUTION IS SPECAC WHAT'S THE PROBLEM?

smiths
A part of Smiths Group plc

Specac

www.specac.com



SPECAC LIMITED,
River House, 97 Cray Ave, St. Mary Cray,
Orpington, Kent BR5 4HE U.K.
T : +44 (0)1689 873134
F : +44 (0)1689 878527
E : sales@specac.com
Registered No. 1008689 England

SPECAC INC.,
50 Sharpe Drive,
Cranston, RI 02920 U.S.A.
T : Toll Free +1 800 447 2558
E : sales@specac.com
W : www.specac.com
Specac is part of Smiths Group plc