

## Specadie

*Produces KBr pellets without using a press*

### Applications

- Halide pellets for IR analysis

### Product Highlights

- No press or pellet holder required
- Rapid sample preparation
- Evacuatable for sample clarity and quality

The evacuatable die produces an 8.5 mm diameter pellet between two highly polished bolt ends.

The die body, which acts as a pellet holder, is mounted on a special Specadie holder



**Specadie Kit 03700**

## Ordering Information Specadie

### Specadie

**03550** Specadie with bolts and seals only

**03700** Specadie Kit  
Includes: Specadie, Specadie Holder, bench mounted wrench, open ended wrench, spare set of bolts, spare set of seals, bottle of KBr powder (50 g)

### Spares and consumables

**03560** Circular Specadie holder

**03570** Set of bolts

**03590** Set of two open ended wrenches, 9/16" (14.3 mm) and 1/2" (12.7 mm) AF

**03595** Bench mounted wrench

**03580** Set of seals

**03610** KBr powder