

<p>Tube, pyrocoated 56UN004 recommended!</p>	9423-393-95091	<p>D6,8 / 5,0 x 28 mm</p> 	56UN001
<p>Tube, uncoated, halfrings</p>	9423-393-95031	<p>D6,8 / 3,6 x 28 mm</p> 	56UN002
<p>Tube, pyrocoated, halfrings</p>	9423-393-95071	<p>D6,8 / 3,6 x 28 mm</p> 	56UN004
<p>Tube, pyrocoated – LONGLIFE thicker coated, recommenden for over night runs</p>	9423-393-95041	<p>D6,8 / 3,6 x 28 mm</p> 	56UN006
<p>Tube, pyrocoated, forked platform recommended for samples with salt matrix</p>	9423-393-95191		<p>56UN007 on demand</p>
<p>Contacts, Zeeman Electrode set of 2 pieces</p>	9423-393-95161	<p>D18,0 / 4,4 x 20,0 mm</p> 	56UN011
<p>Contacts, standard Electrode set of 2 pieces</p>	9423-393-95011	<p>D20 / 5 x 24,4 mm</p> 	56UN012