









Tube, pyrocoated, with rings recommended standard tube, high quality	63-100012-00		56VA002
Tube, pyro coated, with rings, Longlife better coating, recommended for over-night-runs	63-100012-EL	D7.0 / 4.8 x 25.0 mm 	56VAEL2
Plateau-tube, coated also for use with bone platform	63-100011-00		56VA003
Bone-Platform Pyro graphite, inserting in Plateau-tube, not recommended!	63-100013-00	Pyro-Platform 	56VA004
Tube, uncoated, with rings	63-100015-00		56VA005
Tube, coated, with Omega-Platform	63-100026-00		56VA010
Plateau-tube, uncoated also for use with bone platform	63-100014-00		56VA033
Electrodes / contacts for SpectrAA set of 2 pieces	63-100016-00		56VA006
Zeeman-Contacts set of 2 pieces	63-100017-00		56VA007
Shroud for SpectrAA 1 piece	63-100018-00		56VA008
Zeeman Shroud 1 piece	63-100019-00		56VA009